



Aqua Creek Products

THE LEADERS IN RECREATION, FITNESS, AND ABILITY!



No Anchors Required

Item# F-12PPL-HD-AT1

The Patriot Lift is the only completely portable aquatic lift that is independently verified to meet ADA guidelines. It has a lifting capacity of 375 pounds. Our fixing kit (included) provides an ADA solution.

This lift can easily be removed from the pool side when not in use. It features a compact, low profile frame that does not intrude into the pool when not being used. The unique sliding mechanism allows the lift to be deployed and returned to the cart without fastening hardware.




Push lift to pool side, set brakes, slide into place, and it's ready to use.

OPTIONS/ACCESSORIES:

| ITEM # | DESCRIPTION |
|------------|--------------------------------|
| F-120PPC | Cover, standard lockable, blue |
| F-120PPC-T | Cover, standard lockable, tan |
| F-019CA | Cycle attachment |
| F-422PLH-2 | Headrest |
| F-105LAR | Pull-out leg rest |
| F-043SCH | Solar charging station |

Available in a variety of accent colors
(blue plastics standard)



U.S. Patent No.
7,249,386 B2

U.S. Patent No.
D507,769 S

U.S. Patent No.
5,790,995

FEATURES:

- 375 lb (170 kg) weight capacity
- Independently verified to ADA standards
- Rechargeable battery with wall mount charger
- Adjustable footrest & lap belt
- Flip-up arm rests
- Water resistant handset
- Fixing kit included
- 5 year structural warranty
- 5 year pro-rated electronics warranty




PROUDLY SOLD BY:


ADA COMPLIANT

The Patriot Lift

PRODUCT SPECIFICATIONS



APPLICABLE POOL PROFILES



To confirm specific application go to www.aquacreek.com for our **Pool Lift Selector**



Become a member of our new **Dealer and Architect** page. Access to CAD blocks, CSI specs and sales materials.