



AMBASSADOR LIFT

US PATENT NO'S: [7,249,386 B2] [D507,769 S]

MEASUREMENT TAKEN WITH THE LIFT AT FULL LOAD



PRODUCT WEIGHT: 150lbs [68kg]
MAX WEIGHT CAP.: 300lbs [136kg]
MIN ANCHOR SETBACK: 22.5"
MAX ANCHOR SETBACK: 38"
MAX WATER DRAFT: *10"



Inspected And Verified To Meet
 Applicable Sections of ADA
 Guidelines Section 1009.2 By

wyle
 laboratories

A Registered ISO 9001
 Laboratory



9889 GARRYMORE LANE
 MISSOULA, MT 59808
 (406) 549-0769
 FAX (406) 549-2602

NOTE:
 *SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

| | | | |
|-----------|------|------------|-----------|
| PRODUCT: | | AMBASSADOR | |
| PART #: | | F-PPXRD-NA | |
| SCALE: | 1:18 | DATE: | 6/14/2019 |
| DRAWN BY: | KR | REVISION: | 03 |